

## HEMOCRAFT QUEENS NELSON KNIFE

**\$25.00**



**SKU:** PAT-AA5508A



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Hemocraft Queens Nelson Knife - a super handy tool for easy eating! This special knife is made of strong stainless steel and has two cool functions in one.

On one end, it's a sharp knife for cutting food smoothly, like buttering bread or slicing meat. But guess what? The other end is a fork with a cleverly curved shape, making it simple to pick up food without needing another fork.

The handle is designed to feel comfy in your hand, so it's easy to hold and use. It's perfect for folks who might find it hard to use regular cutlery, like those with hand issues or older adults.

With the Hemocraft Queens Nelson Knife, everyone can enjoy their meals with less fuss. It's a practical and empowering kitchen tool that's here to make eating easier and more enjoyable for everyone!

### FEATURES

- The knife can be used with a rocker action involving little movement of the wrist.
- The pronged end can be used as a fork.
- The part of the blade near the handle is not sharpened.
- The stainless steel knife blade and plastic handle.
- Dishwasher safe