

HEMCRAFT RAISED TOILET SEAT

\$189.00



SKU: PAT-091558873



PRODUCT DESCRIPTION

RAISED TOILET SEAT WITH ARMRESTS PURPOSE AND BENEFITS

The Homecraft Raised Toilet Seat with Armrests is Ideal if you Impaired Movement and are looking for a Toilet Aid to Help Increase Your Independence when using the Toilet.

Not only will this Raised Toilet Seat fit over most Standard Toilets, but you can Easily Adjust the Arm Heights to Three Different Heights ranging from 60mm, 100mm and 150mm to Suit Your Needs.

FEATURES

- Features removable, rotating armrests and is height adjustable.
- Easy to clean.
- Fits most standard toilets.
- Secured at the rear of the seat using the existing holes on the toilet basin.
- Replaces existing seat and lid.
- Plastic moulded armrests.

SPECIFICATIONS

- Overall Width: 500mm
- Overall Depth: 475mm

- Overall Height: 380mm
- Seat Height: 60mm, 100mm and 150mm
- Seat Width: 360mm
- Seat Depth: 425mm
- Aperture Size: 210×275mm
- Weight: 4.3kg
- SWL: 150kg

Also have a look at [savanah toilet seat raiser](#).