



HEMOCRAFT ROPE LADDER BED HOIST



SKU: PAT-AA3660

PRODUCT DESCRIPTION

Purpose & Benefits of Hemocraft Rope Ladder Bed Hoist

The Hemocraft Rope Ladder Bed Hoist is designed to assist individuals in safely and independently sitting up in bed. This device is particularly useful for people with limited mobility, such as those recovering from surgery, the elderly, or individuals with disabilities. It provides a means of increasing independence in the bedroom, reducing the need for caregiver assistance for simple positioning.

Features

- **Enhanced Mobility:** Allows users to pull themselves up to a sitting position with minimal assistance, promoting independence and mobility.
- **Simple Installation:** The bottom ends of the cords tie around the bed legs, making it easy to install on most bed types.
- **Durable Construction:** Features robust plastic ladder rungs that are easy to grasp, providing a reliable and sturdy aid for moving to a seated position.
- **Suitable for Various Bed Types:** Compatible with single or double beds, offering versatility in use across different sleeping arrangements.
- **Portable:** Lightweight and easy to remove or transfer to another bed as required, enhancing its utility for home or travel.

Specifications

- **Rung Width:** 195 mm
- **Rung Diameter:** 25 mm
- **Cord Length (to first rung):** 3m

The Homecraft Rope Ladder Bed Hoist is a practical and effective solution for those who require a little extra help sitting up in bed. It's designed to foster self-reliance and improve the quality of life for individuals with reduced mobility, making daily life safer and more accessible.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au