



## HEMOCRAFT SHERWOOD BARIATRIC COMMODE

**\$626.00**



**SKU:** PAT-091422963

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Homecraft Sherwood Bariatric Commode** is designed to deliver stability, comfort, and confidence for larger-bodied users in the home. Featuring a sturdy steel frame with non-slip feet, a thickly padded seat and backrest, and height-adjustable legs, it supports extended use and assists with safe transfers. The removable armrests, hygienic pan and seat cover, and wide aperture promote independence and ease of use.

### Features

- Height-adjustable frame (480–635 mm) for customised comfort and optimal positioning
- Extra-wide, thickly padded seat and backrest (seat width: 660 mm; seat depth: 500 mm)
- Removable armrests facilitate side transfers from a bed or wheelchair
- Comes with a square commode pan and separate seat cover for hygiene and convenience
- Built for bariatric use with a strong steel frame and a maximum user weight of 254 kg

### Specifications

- **Seat Width:** 660 mm

- **Seat Depth:** 500 mm
- **Seat Height Range:** 480-635 mm
- **Aperture Size:** 280 × 230 mm
- **Footprint (Width × Depth):** 670 × 590 mm
- **Armrest Height Above Seat:** 220 mm
- **Maximum User Weight:** 254 kg :