



HEMOCRAFT SHERWOOD PLUS BARIATRIC PERCHING AND SHOWER STOOL



SKU: PAT-091156678-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS OF THE HEMOCRAFT SHERWOOD PLUS BARIATRIC PERCHING AND SHOWER STOOL

The Hemocraft Sherwood Plus Bariatric Perching and Shower Stool is designed to provide sturdy and comfortable seating support for individuals requiring assistance during perching or showering activities. Its robust construction and adjustable features cater to bariatric users, ensuring safety and convenience in daily routines. Ideal for use in bathrooms or other areas where perching support is needed, this stool enhances user stability and comfort.

FEATURES

- **Height-Adjustable Legs:** Equipped with pin clips, the legs can be adjusted to accommodate different user heights, promoting optimal comfort and posture.
- **Non-Slip Feet:** The stool features non-slip feet to prevent movement on slippery surfaces, enhancing user safety during use.
- **Generously Padded Seat:** The extra-wide seat is generously padded and covered with vinyl for additional comfort and easy cleaning.
- **Cross-Braced Construction:** Built with 2.5 cm epoxy-coated steel tubing, the stool provides ultra-strong seating support for bariatric users.
- **Optional Depth-Adjustable Backrest:** The shower stool variant includes a depth-adjustable backrest for added support and comfort.
- **Tubular Armrests:** Features tubular armrests that support an upright posture and assist with transfers.

SPECIFICATIONS

Dimensions	Shower Stool	Perching Stool
Seat width	61cm	52.1cm
Seat depth	43cm	32.4cm
Height adjustment	50.8 - 68.6cm	55.9 - 73.7cm
Width between arms	67cm	61cm
Footprint at highest setting	-	77 x 66cm
Maximum User Weight	255kg	255kg

The Homecraft Sherwood Plus Bariatric Perching and Shower Stool offers a reliable and comfortable solution for individuals seeking sturdy support during perching or showering activities, with features tailored to enhance safety and user experience.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au