

HOMECRAFT SHERWOOD PLUS BARIATRIC PERCHING AND SHOWER STOOL

\$535.00



SKU: PAT-091156678



PRODUCT DESCRIPTION

Purpose & Benefits of the Homecraft Sherwood Plus Bariatric Perching and Shower Stool

The Homecraft Sherwood Plus Bariatric Perching and Shower Stool is designed to provide sturdy and comfortable seating support for individuals requiring assistance during perching or showering activities. Its robust construction and adjustable features cater to bariatric users, ensuring safety and convenience in daily routines. Ideal for use in bathrooms or other areas where perching support is needed, this stool enhances user stability and comfort.

Features

- **Height-Adjustable Legs:** Equipped with pin clips, the legs can be adjusted to accommodate different user heights, promoting optimal comfort and posture.
- **Non-Slip Feet:** The stool features non-slip feet to prevent movement on slippery surfaces, enhancing user safety during use.
- **Generously Padded Seat:** The extra-wide seat is generously padded and covered with vinyl for additional comfort and easy cleaning.
- **Cross-Braced Construction:** Built with 2.5 cm epoxy-coated steel tubing, the stool provides ultra-strong seating support for bariatric users.
- **Optional Depth-Adjustable Backrest:** The shower stool variant includes a depth-adjustable backrest for added support and comfort.
- **Tubular Armrests:** Features tubular armrests that support an upright posture and assist with transfers.

Specifications

Dimensions	Shower Stool	Perching Stool
------------	--------------	----------------

Seat width	61cm	52.1cm
Seat depth	43cm	32.4cm
Height adjustment	50.8 - 68.6cm	55.9 - 73.7cm
Width between arms	67cm	61cm
Footprint at highest setting	-	77 x 66cm
Maximum User Weight	255kg	255kg

The Homecraft Sherwood Plus Bariatric Perching and Shower Stool offers a reliable and comfortable solution for individuals seeking sturdy support during perching or showering activities, with features tailored to enhance safety and user experience.