



**HEMECRAFT SHERWOOD PLUS BARIATRIC
SHOWER STOOL WITH ARMS - PADDED SEAT &
BACK**

\$638.00



SKU: PAT-091156694

PRODUCT DESCRIPTION

Purpose & Benefits

The **Homecraft Sherwood Plus Bariatric Shower Stool** provides a stable, extra-wide seating option for use in the shower. Its durable construction and padded seat offer added comfort and safety for bariatric users.

Features

- Extra-wide design for increased comfort and support
- Generously padded vinyl-covered seat for added comfort
- Cross-braced frame and 2.5cm epoxy-coated steel tubing for enhanced strength
- Height adjustable with pin clips to suit user preference
- Non-slip feet help prevent movement on wet surfaces
- Available with a depth-adjustable backrest for improved support

Specifications

- **Seat Width:** 61 cm

- **Seat Depth:** 43 cm
- **Height Adjustment:** 50.8 – 68.6 cm
- **Width Between Arms:** 67 cm
- **Maximum User Weight:** 255 kg