

## TERRY CLOTH SOCK AID

**\$40.00**



**SKU: PAT-A7273**



## PRODUCT DESCRIPTION

## PURPOSE AND BENEFITS

The Homecraft Terry Cloth Sock Aid makes sock application easy. The flexible, plastic sock aid has slits in the plastic core allowing it to bend around the heel. The inside is lined with nylon to reduce friction. The outside is lined with terry cloth to hold the sock firmly until it is on the foot. Two pulling tabs with loop handles permit one-handed application. Hand wash only.

## FEATURES

- Length 230mm
- Weight 73g