

## HUMAN CARE HEALTHCARE BEDS



**SKU:** 100781

### PRODUCT DESCRIPTION

### VARIATIONS

Available in a standard and plus sized version. Choose between white or wooden boards.

### PURPOSE AND BENEFITS

The FloorLine LTC with its home environment design is the optimal choice for long-term care facilities. Its one-piece design makes it easy to fold and store away efficiently. The FloorLine LTC features all the standard required positions often found in hospital beds such as Trendelenburg, reverse Trendelenburg and cardiac positioning.

The FloorLine LTC is the perfect solution for preventing fall injuries. Its minimum height from the floor is just 9.5 cm(3.7") providing a safe and considerate way of taking care of people who may be at risk of falling out of their bed.