

## IC280 HOMECARE BED



**SKU:** IC280

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Essential with high lift.

The IC280 Bed combines the essential features of a homecare bed, with ultra high lift function enabling easier patient care. The three function design is compatible with all Icare accessories and attachments.

On the standard casters height adjustment: Min – 390mm Max – 810mm

With the larger casters it is 30mm different: Min – 420mm Max – 840mm

Please click [here](#) to download PDF Brochure.

### Available Sizes:

#### (1) LONG SINGLE

Mattress Surface Size: 2030 x 900

Base Size: 2150 x 960

#### (2) KING SINGLE

Mattress Surface Size: 2030 x 1070

Base Size: 2150 x 1120

### **Accessories Available:**

- Four Rail Low Side Safety Rail
- Four Rail High Side Grab Rail
- Bed Stick with Safety Return
- Full Length Bed Rail
- Over Bed Pole
- Over Bed Table

### **FEATURES**

- 180kg Safe Working Load
- Available With Headboards & Footboards
- Rechargeable Battery Backup
- Heavy Duty Castors
- Strong Breathable Mattress Platform
- End And Corner Mattress Retainer Bars
- 5 Year Limited Warranty
- Back Lit Wipeable Hand Piece
- Vertical Lift In Place

### **SPECIFICATIONS**

- Safe working load – 180 kg on long single and king single sizes. The safe working load includes any model of Icare mattresses.
- Hand Piece – Backlit for nighttime use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.
- Fabric – Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Mattress Surface – Strong metal mesh under mattress allowing increased airflow
- No pull motors for added safety.
- ctuators – 1x 6000 newtons – 24 watt head adjust actuator. 1x 6000 newtons – 24 watt knee adjust actuator. 2x 6000 newtons – 24 watt height adjust actuators (2x 8000 newton on double and queen sizes)
- No pull motors for added safety.
- Wheels – 4 heavy duty swivel castors with two lockable castors.
- Head lift angle 62° and knee break angle 42°

### **Colours Available**

[ICare-Upholstery-Colour-Choices.pdf](#)