



## ICARE ACTIVEFLOAT TOPPER OVERLAY

**\$610.00 - \$1,003.00** Price range: \$610.00 through \$1,003.00



**SKU:** ICOF-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Icare ActiveFloat Overlay, utilising cutting-edge technology, employs a highly resilient and exceptionally breathable substance constructed to mitigate shear forces and provide pressure alleviation. This top-grade overlay minimises shear in all directions, enhancing the user's comfort. Shear reduction is achieved through the implementation of SlideSkin technology, which employs an additional membrane allowing the cover to correspond with the movements of the user.

### FEATURES

- **Ventilated Side Fabric** - The side wall of this overlay cover utilises Vant fabric, a three-dimensional material made of polypropylene that is highly breathable and resistant to bacteria, thereby significantly improving air circulation
- **SlideSkin™** technology reduces shear with a second membrane, enabling the cover to move or slide with the users movements
- Suitable for use on some alternating air mattresses. ActiveX™ Overlays are often used to increase comfort and reduce heat loss when used with alternating air mattresses

### SPECIFICATIONS

- **Long Single:** 2030mm x 900mm x 50mm - Min. 30kg | Max. 215kg
- **King Single:** 2030mm x 1070mm x 50mm - Min. 30kg | Max. 215kg
- **Short Double:** 1900mm x 1350mm x 50mm - Min. 30kg | Max. 215kg
- **Queen:** 2030mm x 1520mm x 50mm - Min. 30kg | Max. 215kg