

## ICARE ACTIVEX MATTRESS BOLSTER

**\$198.00 - \$288.00** Price range: \$198.00 through \$288.00



**SKU:** MB-H

## PRODUCT DESCRIPTION

Suited for IC333EL Bed.

Two choices of feel – IC15 Firm or IC20 Medium Mattress.

Bolster size is 300mm and is compatible for IC333EL

How pressure care works: The ActiveX Mattresses are heat and pressure sensitive which means that as it gets warmer it becomes softer. So when there is a pressure point developing which usually gets hot and enflamed, this area of the mattress will become softer which allows the blood to flow freely again.