

ICARE ACTIVEX MATTRESS BOLSTER

\$178.00 - \$259.00Price range: \$178.00 through \$259.00





SKU: MB-H

PRODUCT DESCRIPTION

Suited for IC333EL Bed.

Two choices of feel - IC15 Firm or IC20 Medium Mattress.

Bolster size is 300mm and is compatible for IC333EL

How pressure care works: The ActiveX Mattresses are heat and pressure sensitive which means that as it gets warmer it becomes softer. So when there is a pressure point developing which usually gets hot and enflamed, this area of the mattress will become softer which allows the blood to flow freely again.