


## ICARE ACTIVEX TOPPER

**\$384.00 - \$544.00** Price range: \$384.00 through \$544.00



**SKU:** ICT-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Enhance pressure relief and comfort on your existing mattress with the innovative Icare ActiveX Pressure Care Overlay. Specially designed, this overlay is heat and pressure sensitive, adapting to the body's contours and providing softness around pressure points and inflamed areas.

Crafted from a high-grade elastic foam known as ActiveX, this overlay responds to the temperature around inflamed areas. As it becomes warmer, ActiveX becomes softer, ensuring a gentle and supportive surface around pressure points. This feature promotes proper blood circulation, allowing for improved comfort and well-being. Experience the benefits of the Icare ActiveX™ Pressure Care Overlay, a solution designed to enhance pressure relief and optimize comfort for a restful sleep experience.

### FEATURES

- **Compatible with Alternating Air Mattresses:** Suitable for use with alternating air mattresses to enhance comfort and reduce heat loss.
- **Pressure Redistribution:** The ActiveX 84 material responds to temperature, softening around pressure points to improve blood circulation.
- **Improved Blood Circulation:** The overlay promotes better circulation by providing a softened surface around pressure points.

### SPECIFICATIONS

#### SWL

Minimum 30kg (66 lbs)  
Maximum 215kg (474 lbs)

**Depth**

5cm (1.9 inches)

**Mattress size**

Short Double 1900 x 1350mm

Long single 2030 x 900mm

Queen 2030 x 1520mm

King single 2030 x 1070mm