

**ICARE BLANKET CRADLE****\$107.00****SKU: ACOBC****PRODUCT DESCRIPTION**

Bed mounted blanket cradle. The Icare Blanket cradles is designed to support the weight of bed linen. This can relieve the lower legs and feet from pressure while also keeping sheets or blankets away from sores or wounds. Blanket cradles can also promote air circulation to help keep skin dry.

Suits all Icare beds (*except Icare Companion base*) Charcoal colour, easy to mount to the bed base/mattress platform, moves with the mattress platform.

Easy height adjustment system, and removable. Measured from mattress platform 50cm to 43cm.

Width 420mm