




ICARE C777 PREMIUM HOMECARE BED



SKU: IC777

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Bring beauty and character into your room with exceptional bed functions, made in Australia. This four function bed, features dual lift motors totalling 12,000 newtons lifting power. It also has unlimited under-bed clearance and a cordless remote with preset functions.

Please [click here](#) to download full brochure.

CODES & SIZING

Colours Available

[ICare-Upholstery-Colour-Choices.pdf](#)

KEY FUNCTIONS

SPECIFICATIONS

- 200kg safe working load on long single and king single sizes and 300kg for double and queen sizes. The safe working load includes any model of ICare mattresses.
- Long-life rechargeable battery back up, lasting 2-5 days of usual use.
- Easy press button remote, backlit when activated for use. Wipeable.
- Cordless remote with flashlight, and preset functions.

- Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Strong metal mesh mattress platform allowing increased airflow
- Actuators – 1x 6000 newton – 24 watt head adjust actuator. 1x 6000 newton – 24 watt knee adjust actuator. 2x 6000 newton – 24 watt height adjust actuators (2x 8000 newton on double and queen sizes)
- No pull motors for added safety.
- Wheels – 8x heavy duty swivel castors with one lockable castor on each corner.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au