




ICARE C777 PREMIUM HOMECARE BED



SKU: IC777

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Bring beauty and character into your room with exceptional bed functions, made in Australia. This four function bed, features dual lift motors totalling 12,000 newtons lifting power. It also has unlimited under-bed clearance and a cordless remote with preset functions.

Please [click here](#) to download full brochure.

CODES & SIZING

Colours Available

[ICare-Upholstery-Colour-Choices.pdf](#)

KEY FUNCTIONS

SPECIFICATIONS

- 200kg safe working load on long single and king single sizes and 300kg for double and queen sizes. The safe working load includes any model of ICare mattresses.
- Long-life rechargeable battery back up, lasting 2-5 days of usual use.
- Easy press button remote, backlit when activated for use. Wipeable.
- Cordless remote with flashlight, and preset functions.

- Hard wearing commercial grade fabric. Vinyl is treated with antibacterial properties.
- Strong metal mesh mattress platform allowing increased airflow
- Actuators – 1x 6000 newton – 24 watt head adjust actuator. 1x 6000 newton – 24 watt knee adjust actuator. 2x 6000 newton – 24 watt height adjust actuators (2x 8000 newton on double and queen sizes)
- No pull motors for added safety.
- Wheels – 8x heavy duty swivel castors with one lockable castor on each corner.