



## ICARE EXTENDABLE BED SIDE RAIL

**\$236.00**



**SKU: ACESR**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Icare Extendable Side Rail is designed to reduce the risk of bed users from falling, or accidentally rolling out. The rail can be retracted allowing users to easily access and exit the bed without obstruction. They are particularly useful in situations where the side rail is not needed at all times.

### FEATURES

- Charcoal colour
- Slip-on attach system with four rails for maximum safety
- Position anywhere along either side of the bed where required
- The extension can be utilized in either direction, towards the head or the feet
- Locks in two positions, standard width or fully extended

### SPECIFICATIONS

- Width 690 mm – 950 mm
- Height above mattress platform 450mm
- Space between rails 100mm
- Max user weight 200kg