



ICARE EXTENDABLE BED SIDE RAIL



SKU: ACESI

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Icare Extendable Side Rail is designed to reduce the risk of bed users from falling, or accidentally rolling out. The rail can be retracted allowing users to easily access and exit the bed without obstruction. They are particularly useful in situations where the side rail is not needed at all times.

FEATURES

- Charcoal colour
- Slip-on attach system with four rails for maximum safety
- Position anywhere along either side of the bed where required
- The extension can be utilized in either direction, towards the head or the feet
- · Locks in two positions, standard width or fully extended

SPECIFICATIONS

- Width 690 mm 950 mm
- · Height above mattress platform 450mm
- Space between rails 100mm
- Max user weight 200kg

For further information go to medimart.com.au or give us a call: 1300 633 063



SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au