


ICARE FULL LENGTH FOLD-DOWN BED RAILS

\$224.00 - \$448.00 Price range: \$224.00 through \$448.00



SKU: ACSFR

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Full length fold-down bed rails, available in left and right side. Your choice of the colour finish allows for a clean finish in all environments. Easy installation yet an incredibly strong structure means you can feel safe and secure without sacrificing space, aesthetics or comfort. **For ICARE Beds ONLY.**

- Compatible with all Icare beds
- Fold down design
- 4 Rails for maximum safety
- Slip on attach system
- Position anywhere along to side as needed
- Charcoal colour
- Maximum user weight: 200kg
- Width: 1385mm
- Height above mattress platform: 495mm
- Space between rails: 80mm