


## **ICARE FULL LENGTH FOLD-DOWN BED RAILS**

\$196.00 - \$392.00Price range: \$196.00 through \$392.00





SKU: ACSFF

## PRODUCT DESCRIPTION

## **PURPOSE AND BENEFITS**

Full length fold-down bed rails, available in left and right side. Your choice of the colour finish allows for a clean finish in all environments. Easy installation yet an incredibly strong structure means you can feel safe and secure without sacrificing space, ascetics or comfort. **For ICARE Beds ONLY.** 

- Compatible with all Icare beds
- Fold down design
- 4 Rails for maximum safety
- Slip on attach system
- Position anywhere along to side as needed
- Charcoal colour
- Maximum user weight: 200kg
- Width: 1385mm
- Height above mattress platform: 495mm
- Space between rails: 80mm