



ICARE FULL LENGTH SIDE RAIL COVER

\$205.00



SKU: ACSRCFL

PRODUCT DESCRIPTION

Features and Benefits

The iCare Low Side Rail Cover is specifically designed to enhance user safety and comfort when using iCare low bed rails. Featuring full-height padding and constructed with soft, medical-grade Zerotec fabric, it offers a protective barrier that reduces the risk of bumps or entrapment.

This hygienic, easy-to-clean padded sleeve slips securely over the rail and fastens with extra-strong Velcro, ensuring a snug, stable fit throughout use. Perfect for users at risk of injury from involuntary movements, especially in aged care and home care environments.

Ideal for improving patient comfort while maintaining high hygiene standards, this accessory is a valuable addition to any iCare bed setup.

Key Features

- Full-height padded cover for iCare low bed rails
- Provides cushioning to minimise impact and injury
- Soft, medical-grade Zerotec fabric
- Wipeable, hygienic surface for easy cleaning
- Secure Velcro closure ensures a stable fit
- Simple to install and remove

- Suitable for aged care and home care environments

Specifications

- **Width:** 690mm
- **Height above mattress platform:** 422mm
- **Material:** Zerotec two-way stretch medical-grade fabric
- **Attachment:** Slip-on with Velcro fastening
- **Sold as:** Pair

Estimated Shipping Dimensions

- **Packed Length:** 450mm
- **Width:** 250mm
- **Height:** 50mm
- **Weight:** Approx. 0.7kg