

**ICARE HYDROSENSE BED PAD****\$330.00****SKU: PPSMP****PRODUCT DESCRIPTION****PURPOSE AND BENEFITS**

The essential cooling element

The mattress pad provides an essential all night temperature regulating, can be used for long periods of time.

Size: 190 x 90cm, designed for use by one person.

This item must be used in conjunction with the [HydroSense Unit](#).

FEATURES

- Cools to 8°C and heats to 48°C
- The machine has night settings and a built in timer, to prevent chills.
- 80 Watt power consumption
- No electro-magnetic fields
- Whisper quiet motor
- HydroSense Bed Pad with silicone medical grade tubing
- Available in one size: Single 190 x 90cm
- A range of accessories available including small and medium size multi use pads, back and shoulder vest, neck, limb and eye wrap, foot and leg cuff, and lumbar strap

[Click here to view the Hydrosense Unit](#)