

## ICARE BEDS - IC333 HOMECARE BED



**SKU:** IC333

## PRODUCT DESCRIPTION

### ICARE IC333 HOMECARE BED BENEFITS

The ICare IC333 is an electrically operated bed that has been designed for both home care and [healthcare bed](#) facilities. This bed has adjustable back rest, knee break, is height adjustable and also offers the user trendelenburg tilt. The IC333 also uses a flexible slat systems that support the mattress and comes with an upholstered surround. Available in Onyx or Stone fabric.

Companion Beds and head boards also available.

### Key Features

- Maximum load rated 200kg
- Upholstered head board to match IC333
- Back rest and knee rest
- Trendelenburg tilt 15 degrees
- Heavy-duty locking castors
- Choice of fabrics for upholstered
- Dual lift motors totalling 12,000 Newtons lifting power

### Size Options

- **Long Single:** 2030mm x 900mm
- **Extra Long Single:** 2330mm x 900mm
- **King Single:** 2030mm x 1070mm
- **Queen:** 2030mm x 1500mm
- **Split Queen:** 2030mm x 740mm per side

- **Split King (parter set-up):** 2030mm x 900mm per side

## Colour & Fabric Options

View the full range of options here - [ICare-Upholstery-Colour-Choices.pdf](#)

## Accessories

- High side rail
- Low side rail
- Bed sticks
- Full length fold down bed rail
- Over bed table
- Over bed pole

## Specifications

- Safe Working Load:
  - 200kg for Long Single & King Single (Includes any Icare mattress)
  - 300kg for Double & Queen (Includes any Icare mattress)
- Height Range:
  - Minimum Height: 22cm (floor to top of mattress platform)
  - Maximum Height: 66cm
- Lift Angles:
  - Head Lift Angle: 62°
  - Knee Break Angle: 42°
- Motors & Actuators:
  - 1x 6000N / 24W Head Adjust Actuator
  - 1x 6000N / 24W Knee Adjust Actuator
  - 2x 6000N / 24W Height Adjust Actuators (2x 8000N for Double & Queen)
  - No-Pull Motors for added safety
- Hand Control Unit:
  - Backlit for Night Use
  - Wired Stretch Cord
  - Hook Attachment for Side of Bed
  - Wipeable for Hygiene
- Castors:
  - 8 Heavy-Duty Swivel Castors
  - One Lockable Castor on Each Corner

Buy the ICare Bed IC333 at Medimart Australia, as we have [NDIS Payment method](#) available.