

ICARE BEDS - IC333 HOMECARE BED



SKU: IC333

PRODUCT DESCRIPTION

ICARE IC333 HOMECARE BED BENEFITS

The ICare IC333 is an electrically operated bed that has been designed for both home care and [healthcare bed](#) facilities. This bed has adjustable back rest, knee break, is height adjustable and also offers the user trendelenburg tilt. The IC333 also uses a flexible slat systems that support the mattress and comes with an upholstered surround. Available in Onyx or Stone fabric.

Companion Beds and head boards also available.

Key Features

- Maximum load rated 200kg
- Upholstered head board to match IC333
- Back rest and knee rest
- Trendelenburg tilt 15 degrees
- Heavy-duty locking castors
- Choice of fabrics for upholstered
- Dual lift motors totalling 12,000 Newtons lifting power

Size Options

- **Long Single:** 2030mm x 900mm
- **Extra Long Single:** 2330mm x 900mm
- **King Single:** 2030mm x 1070mm
- **Queen:** 2030mm x 1500mm
- **Split Queen:** 2030mm x 740mm per side

- **Split King (parter set-up):** 2030mm x 900mm per side

Colour & Fabric Options

View the full range of options here - [ICare-Upholstery-Colour-Choices.pdf](#)

Accessories

- High side rail
- Low side rail
- Bed sticks
- Full length fold down bed rail
- Over bed table
- Over bed pole

Specifications

- Safe Working Load:
 - 200kg for Long Single & King Single (Includes any Icare mattress)
 - 300kg for Double & Queen (Includes any Icare mattress)
- Height Range:
 - Minimum Height: 22cm (floor to top of mattress platform)
 - Maximum Height: 66cm
- Lift Angles:
 - Head Lift Angle: 62°
 - Knee Break Angle: 42°
- Motors & Actuators:
 - 1x 6000N / 24W Head Adjust Actuator
 - 1x 6000N / 24W Knee Adjust Actuator
 - 2x 6000N / 24W Height Adjust Actuators (2x 8000N for Double & Queen)
 - No-Pull Motors for added safety
- Hand Control Unit:
 - Backlit for Night Use
 - Wired Stretch Cord
 - Hook Attachment for Side of Bed
 - Wipeable for Hygiene
- Castors:
 - 8 Heavy-Duty Swivel Castors
 - One Lockable Castor on Each Corner

Colour & Fabric Options

- Premium Fabric: Bison (Vinyl)
- Premium Fabric: Mist

- Premium Fabric: Noir
- Premium Fabric: Royal (Discontinued)
- Premium Fabric: Sand (Vinyl)
- Premium Fabric: Slate (Vinyl)
- Onyx
- Silver
- Smoke
- Stone
- Vinyl Dove
- Vinyl Ebony

Buy the ICare Bed IC333 at Medimart Australia, as we have [NDIS Payment method](#) available.

For further information go to **[medimart.com.au](https://www.medimart.com.au)** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au