



## ICARE IC25 ACTIVEX SOFT MATTRESS

**\$750.00 - \$1,419.00** Price range: \$750.00 through \$1,419.00



**SKU:** ICX-H

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

For clients who like a soft mattress, or who spend a long time in bed. The feel of the IC25 Mattress is best suited to aged or petite patients.

The IC25 Mattresses have almost no partner disturbance or roll together.

### KEY FUNCTIONS

- **Pressure Redistribution**

The IC25 mattress features the ActiveX 84 material is the main comfort layer in this mattress. ActiveX is a high grade of elastic type foam which is responsive to temperature around inflamed areas. ActiveX will become softer as it becomes warmer. So around a pressure point it will always soften and allow the blood to circulate freely.

- **Solid Support Foam Base**

XP Support foam is a high-density foam that is used in the bottom layer of the IC25. XP Support foam is an open cell material that breaths to allow air circulation but is high density to prevent bottoming out. This material will easily contour to any profiling bed.

- **Temperature Regulating Cover**

The IC25 cover is made with the new Bluetec cooling fabric. This is all way stretch fabric to conform with body and mattress contour. It is under woven with a cooling yarn to regulate body temperature.

- **Ventilated Base Fabric**

Vant fabric is a three-dimensional polypropylene material that has high air permeability and is anti-bacterial. We use this on the side wall of the IC25 mattress covers to greatly increase air movement.

- **Medical Covers Available**

Zerotec fabric is a medical barrier fabric that's has high elasticity in every direction which reduces shear on pressure redistribution surfaces. It is waterproof and wipeable but has microporous breathable properties. This cover has welded seams and waterfall flaps for fluid control.

## SPECIFICATIONS

- **Safe Working Load:** 30kg - 200kg

	Height	Length	Width
Long Single:	20cm	203cm	90cm
King Single:	20cm	203cm	107cm
Long Double:	20cm	203cm	135cm
Queen:	20cm	203cm	153cm
Split Queen:	20cm	203cm	75cm