

ICARE IC280 JUNIOR HOMECARE BED



SKU: IC280JNR



PRODUCT DESCRIPTION

Purpose & Benefits of the ICare IC280 Junior Homecare Bed

The ICare IC280 Junior Homecare Bed is designed specifically for younger users, offering a smaller frame and enhanced comfort for children or individuals requiring a compact-sized bed in a homecare environment. With electric height and positioning adjustments, this bed provides exceptional support for caregivers, allowing for easier patient care while ensuring the comfort and safety of the user. The IC280 Junior is perfect for homecare settings, offering a blend of functionality and style to suit various needs.

Features

- **150 kg Safe Working Load:** The bed is designed to support users up to 150 kg, ensuring safety and durability during long-term use.
- **Strong Breathable Mattress Platform:** The breathable mattress platform promotes airflow, enhancing user comfort and extending mattress longevity.
- **Optional Headboard and Footboard:** Customise the bed with an optional headboard and footboard to suit individual needs and room aesthetics.
- **End and Corner Mattress Retainer Bars:** These retainers help keep the mattress securely in place, providing additional safety and preventing unwanted movement.
- **5-Year Warranty:** Comes with a 5-year warranty, offering peace of mind and assurance of long-term reliability.
- **External Rechargeable Battery Back-up:** In case of a power outage, the external battery back-up ensures that essential bed functions remain operational.
- **Backlit Wipeable Hand Piece:** The bed's hand control features a backlit design for easy use in low-light conditions, and its wipeable surface allows for quick and hygienic cleaning.

Specifications

- **Bed Size:** 1800 mm x 860 mm
- **Colours:** Blush; Sky; Kiwi; Scarlett;
- **Weight Capacity:** 150 kg
- **Max Height:** 78 cm from floor to top of mattress platform
- **Minimum Height:** 37 cm from floor to top mattress platform
- **Head list angle:** 62 degrees
- **Knee Break Angle:** 42 degrees

Colours Available

[ICare-Upholstery-Colour-Choices.pdf](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au