


## ICARE IC30 JUNIOR GEL MATTRESS



**SKU:** IC30JNR

## PRODUCT DESCRIPTION

### Purpose & Benefits of the ICare IC30 Paediatric Gel Mattress

The ICare IC30 Paediatric Gel Mattress is specifically designed to provide exceptional comfort, support, and pressure relief for children in home care settings. Incorporating gel-infused ActiveX™ material, this mattress helps regulate temperature, ensuring a cool and comfortable sleep surface. The addition of charcoal aids in absorbing odours and moisture, promoting a clean and hygienic sleeping environment. Its medium firmness offers an optimal balance between support and comfort, making it ideal for paediatric users requiring pressure care.

### Features

- **Gel-Infused ActiveX™ Technology:** Provides temperature regulation and pressure relief, enhancing comfort and support for paediatric users.
- **Charcoal-Infused Material:** Absorbs odours and moisture, maintaining a hygienic sleep environment.
- **Firm Edge Support:** Maintains the mattress's shape and structure, reducing the risk of roll-out and ensuring safety.
- **Medium Firmness:** Offers a balanced feel, providing both comfort and necessary support for growing children.
- **Bluetec Cooling Cover:** Features a temperature-regulating cover that enhances the cooling effect, ensuring a comfortable sleep surface.

### Specifications

- **Size:** 1800 mm x 800 mm
- **Mattress Height:** 200 mm

- **Weight Capacity:** Suitable for users up to 200 kg
- **Risk Category:** High Risk
- **Warranty:** 10-Year Limited Warranty
- **Cleaning:** Bluetec cover is machine washable for easy maintenance
- **Safe Working Load:** 30kg - 200kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)