



ICARE IC333EL HOMECARE BED



SKU: IC333EL

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Icare IC333EL is an elegant and functional homecare bed.

The IC333EL Extra Long Bed has the same functions, features and durability as the standard size IC333, just longer. The extra length is added at the foot end.

Mattress Bolster are finished in the same fabric as the ActiveX™ Mattresses.

Bolsters available in: IC15 or IC20 Mattress

Overall Bed Length: 2450mm Extension Bolster Length: 300mm

FEATURES

- 200kg Safe Working Load On Long Single And King Single
- Strong Breathable Mattress Platform
- Lockable Heavy Duty Castors
- Available With Headboards
- External Rechargeable Battery Backup
- 300kg Safe Working Load On Double And Queen

- End And Side Mattress Retainer Bars
- 10 Year Limited Warranty
- Back Lit Wipeable Hand Piece
- No Gap Between Mattresses On King Beds (2x Long Singles)
- Unlimited Under Bed Clearance
- Ultra Low Design Reducing Risk Of Falls

SPECIFICATIONS

	Long Single	King Single	Double	Queen
Overall Height:	22 to 66cms	22 to 66cms	22 to 66cms	22 to 66cms
Overall Length:	245cm	245cm	245cm	245cm
Overall Width:	96cm	112cm	140cm	158cm
Mattress Surface:	233 x 90cms	233 x 107cms	233 x 135cms	233 x 150cms
Head Lift Angle:	62°	62°	62°	62°
Knee Break Angle:	42°	42°	42°	42°
Weight Capacity:	200kg	200kg	300kg	300kg

Colours Available

[ICare-Upholstery-Colour-Choices.pdf](#)