



ICARE IC555 BARIATRIC BED BASE - KING SINGLE

\$4,621.00 - \$4,784.00 Price range: \$4,621.00 through \$4,784.00



SKU: IC555KS SILVER-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **ICare IC555 Bariatric Bed** is a premium four-function electric bed designed to provide comfort, safety, and durability for bariatric users. Engineered with a high weight capacity and advanced adjustability, it enhances user independence and caregiver efficiency in both home and care facility settings. With a modern design and quiet motor operation, the IC555 offers exceptional stability and support for long-term care needs.

Features

- Four-function electric bed with hi-lo, backrest, knee break, and Trendelenburg/reverse Trendelenburg adjustments
- High weight capacity suitable for bariatric users, ensuring reliable performance and stability
- Quiet, smooth motor operation for comfortable and discreet use
- Heavy-duty steel frame construction for superior strength and durability
- Compatible with ICare headboards, side rails, and accessories for customisation
- Designed for easy assembly and transport with a stylish, low-profile appearance

Specifications

Specification	Details
Model	IC555 Bariatric Bed
Safe Working Load	300 kg
Mattress Platform Size	1070 mm x 2030 mm
Height Range	370 mm – 800 mm
Backrest Adjustment	70°
Knee Break Adjustment	30°
Trendelenburg / Reverse Trendelenburg	12° / 12°
Frame Construction	Reinforced powder-coated steel
Warranty	10 years on frame, 2 years on motors