


ICARE IC555 BARIATRIC BED



SKU: IC555

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The IC555 Model is a premium Bariatric Bed with a heavy weight rating and extra features including vascular support and the option of vinyl or fabric.

FEATURES

- Features head lift and knee break
- Strong metal mesh under mattress allowing increased airflow
- Ultra low design reducing risk of falls
- Vascular support
- Lockable heavy duty castors
- End and side mattress retainer bars
- Backlit wipeable hand piece
- Rechargeable battery backup
- Compatible with side rails and accessories
- Upholstery: Fabric – hard wearing commercial grade | Vinyl – treated with antibacterial properties and waterproof
- Wheels – 8x heavy duty swivel castors with one lockable castor on each corner.
- Head lift angle 62°
- Knee break angle 42°
- Max height – 66cm from floor to top of base.
- Min height – 32cm from floor to top of base.

SPECIFICATIONS

Size	King Single
Mattress Surface	2030mm x 1070mm
Base	2150mm x 1120mm
Height	320mm - 660mm
SWL	300kg

Colours Available

[ICare-Upholstery-Colour-Choices.pdf](#)