


ICARE IC555 BARIATRIC BED



SKU: IC555

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The IC555 Model is a premium Bariatric Bed with a heavy weight rating and extra features including vascular support and the option of vinyl or fabric.

FEATURES

- Features head lift and knee break
- Strong metal mesh under mattress allowing increased airflow
- Ultra low design reducing risk of falls
- Vascular support
- Lockable heavy duty castors
- End and side mattress retainer bars
- Backlit wipeable hand piece
- Rechargeable battery backup
- Compatible with side rails and accessories
- Upholstery: Fabric – hard wearing commercial grade | Vinyl – treated with antibacterial properties and waterproof
- Wheels – 8x heavy duty swivel castors with one lockable castor on each corner.
- Head lift angle 62°
- Knee break angle 42°
- Max height – 66cm from floor to top of base.
- Min height – 32cm from floor to top of base.

SPECIFICATIONS

Size	King Single
Mattress Surface	2030mm x 1070mm
Base	2150mm x 1120mm
Height	320mm - 660mm
SWL	300kg

Colours Available

[ICare-Upholstery-Colour-Choices.pdf](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au

