



ICARE M3 MEDICAL MATTRESS



SKU: ICM3-H

PRODUCT DESCRIPTION

ICare M3 Medical Mattress suited for high risk with ActiveFloat.

The M3 has an unparalleled user comfort rating which is a crucial factor and often overlooked when choosing a suitable mattress.

- **ActiveFloat Layer**

Our latest technology is an extremely elastic and highly breathable material designed to move in every direction even after full compression. This not only reduces shear in every direction but greatly increases user comfort.

- **Edge to Edge**

The full ActiveFloat layer provides edge to edge support for users eliminating ridges and increasing pressure redistribution area.

- **Three Layer Construction**

One piece U shape firm XP Support Foam base and side wall support with medium feel lateral castellated foam inlay. Designed especially to reduce roll-out risk and bed falls.

- **ActiveFloat Layer**

Our latest technology is an extremely elastic and highly breathable material designed to move in every direction even after full compression. This not only reduces shear in every direction but increases user comfort.

- **Zerotec Cover**

Zerotec fabric is a medical barrier fabric that has high elasticity in every direction which reduces shear on pressure redistribution surfaces. It is waterproof and wipeable but has microporous breathable properties. This cover has welded seams and waterfall flaps for fluid control.

- **SlideSkin Anti-Shear**

The unique composition of SlideSkin technology reduces shear with a second membrane. This

enables the cover to move or slide with the users movements.

SPECIFICATIONS

- **Safe Working Load:** 30kg - 215kg

	Height	Length	Width
Long Single:	20cm	203cm	90cm
King Single:	20cm	203cm	107cm
Long Double:	20cm	203cm	135cm
Queen:	20cm	203cm	153cm
Split Queen:	20cm	203cm	75cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au