

ICARE MEDICAL CUSHION COVER

\$47.00



SKU: ICC455



PRODUCT DESCRIPTION

Purpose & Benefits

The Icare Medical Cushion Cover is a high-quality, washable protective cover designed to fit standard 45 × 45 × 7 cm wheelchair or support cushions. Made from breathable, low-shear medical-grade fabric, it helps protect cushions from moisture, bacteria, and everyday wear while maintaining skin integrity and comfort.

Ideal for clinical, aged care, or home environments, the cover features a water-resistant surface and a zip closure for a secure fit. It's easy to remove and machine wash, making it a practical solution for maintaining hygiene and extending cushion life.

Features

- Medical-grade, breathable two-way stretch fabric
- Low-shear surface to reduce friction and skin breakdown
- Water-resistant and machine washable
- Zip closure for secure and easy application
- Suitable for aged care, clinical, or home use

- Designed to fit standard size cushions

Specifications

- Dimensions: 45 x 45 x 7 cm
- Material: Breathable, two-way stretch fabric
- Closure: Zipper
- Water-Resistant: Yes
- Washable: Machine washable
- Compatibility: Fits Icare cushions and similar-sized products