



## ICARE OVER BED POLE

**\$282.00**



**SKU: ACPUB**

## PRODUCT DESCRIPTION

Freestanding Pull up Bar. The Icare Over Bed Pole helps with turning, patient transfer, and bed care. A stable design with heavy-duty steel, compatible with all Icare Beds.

Height from floor to the top of the rail - 190cm Floor Foot Print: H design > Width under bed 75cm > Length under bed 90cm Rail at the top from back of stand to tip of rail - 71cm.

Maximum weight applied 100kg.