


ICARE THERMOPAD BED

\$456.00



SKU: ICTHRM190



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Reduce core body temperature.

The Icare ThermoPad Overlay is designed to reduce core body temperature. Made up of gel infused ActiveX with an embedded cooling gel pad. The Bluetec cover is under woven with a cooling yarn to assist with body temperature. Suitable for use on an existing mattress or with profiling beds.

SPECIFICATIONS

- One Size 90 x 190cm