

ICARE VERSAFLO3 MATTRESS SYSTEM

\$2,423.00 - \$2,656.00 Price range: \$2,423.00 through \$2,656.00



SKU: ICVF3

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The VersaFlo3 has been designed to combine high pressure care needs and comfort. Using a three layer technology, configured to create a comfortable sleep support surface while providing high pressure relief. The mattress foundation is an XP Support layer that reduces over immersion while enabling the mattress to profile. The mattress core uses adjustable cylinders that alternate air pressure to adjust the surface pressure points. ActiveX is used as the primary comfort layer to increase comfort and disperse pressure.

FEATURES

- Three layer technology
- Alternating or Static air function
- Zerotec mattress cover
- Adjustable air pressure
- High risk pressure care
- CPR release tag

Please [Click Here](#) to download VersaFlo3 brochure.

Please [Click Here](#) to download VersaFlo Pressure Care Information brochure.

SPECIFICATIONS

Long Single

- Mattress Size: 2030 x 900

King Single

- Mattress Size: 2030 x 1070

Versaflo 3 Mattress System - Air Pump

- Circulation Intervals: 5min, 10min, 15min, 20min, static
- Air Capacity: 1.5~8L
- Air Pressure: 300-1200H²O
- User Weight: 20~140kg

Versaflo 3 Mattress System - Mattress

- Material: Polyurethane, TPU + Fabric
- Cover: Waterproof
- Air Pressure: Level 4~5
- Safe Working Load: 150kg