



## IKON E TRAVEL FOLDING - GREEN



**SKU:** IK4QPG

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Lightweight and travel-friendly, this mobility scooter folds easily for transport and storage. It features a height-adjustable T-handlebar, comfortable contoured seat with fold-down backrest, and intuitive finger-operated controls. Battery-powered forward/reverse drive, variable speed, and rear-wheel drive ensure smooth performance, while the rechargeable lithium-ion battery with Tru-View gauge and electromagnetic brake provide reliable power and safety. Perfect for home use or taking on an aeroplane.

### FEATURES

- Folds up for transport and storage
- Height adjustable pull handle
- Safe to take on an aeroplane
- Comfortable contoured seat
- Fold down backrest
- Adjustable height T handlebar
- Finger operated controls
- Battery powered indicator
- Forward / Reverse drive settings
- Variable speed setting
- Rear-wheel drive
- Electromagnetic brake
- Rechargeable Lithium-ion battery
- Graduated tru-view Lithium battery gauge

**\*Please note: the bag is sold separately. [Click here to view the Ikon eTravel Transport Bag](#)**

## SPECIFICATIONS

| Specification                   | Details                                         |
|---------------------------------|-------------------------------------------------|
| Overall Length                  | 910 mm                                          |
| Overall Height                  | 870 mm                                          |
| Overall Width                   | 420 mm                                          |
| Folded Length                   | 710 mm                                          |
| Folded Width                    | 420 mm                                          |
| Folded Height                   | 280 mm                                          |
| Turning Circle                  | 1420 mm                                         |
| Seat Width                      | 330 mm                                          |
| Seat Depth                      | 250 mm                                          |
| Seat Height                     | 390 mm (from footrest)                          |
| Pull Handle Height              | 690 - 1000 mm                                   |
| Tiller Height                   | 700 - 750 mm                                    |
| Speed                           | 6.4 km/h                                        |
| Maximum Range                   | Up to 17 km                                     |
| Maximum Climb Angle             | 6°                                              |
| Motor                           | 180 W high-torque motor                         |
| Battery                         | 24V / 9.9Ah Lithium-Ion (supplied)              |
| Battery Charger                 | 2 Amp Lithium-Ion (supplied)                    |
| Tyres                           | Flat-free tyres                                 |
| Front Tyres                     | Ø165 mm                                         |
| Rear Tyres                      | Ø185 mm                                         |
| Product Weight                  | 19 kg (with battery)<br>17 kg (without battery) |
| Maximum User Weight             | 115 kg (on flat ground)                         |
| Warranty - Structural Frame     | 5 years                                         |
| Warranty - Drive Train Assembly | 18 months                                       |
| Warranty - Other Parts          | 12 months                                       |
| Warranty - Lithium Battery      | 18 months                                       |
| Warranty - Charger              | 12 months                                       |

Color: Green

[Brochure](#)

## OPTIONAL ACCESSORY

[Ikon eTravel Transport Bag](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)