



IKON ENDURO POWER LITE SIDE FOLDING CARBON FIBRE ELECTRIC WHEELCHAIR

\$4,429.00 Original price was: \$4,429.00.\$3,580.00 Current price is: \$3,580.00. Sale Available Online Only!



SKU: IK2PCS-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **ENDURO PowerLite Side Folding Electric Wheelchair** combines durability and convenience with its carbon fibre frame and foldable design for easy transport and storage. It features a comfortable padded seat and backrest, with a fold-down backrest that includes a handy storage pocket. Lift-up armrests and an onboard controller with an LCD display and fingertip controls ensure ease of use. The wheelchair boasts a battery power indicator, variable speed settings, and rear-wheel drive for smooth operation.

Enjoy a comfortable ride with front-wheel suspension, an electromagnetic brake, and anti-tip wheels. The airplane-safe, rechargeable Lithium-ion battery provides reliable power. Additional features include a calf strap, lap belt, and horn for added safety and convenience.

FEATURES

- Carbon fibre frame
- Folds up for transport and storage
- Comfortable padded seat and backrest
- Fold down backrest with storage pocket
- Lift up armrests
- Onboard controller with LCD display and fingertip controls Battery power indicator
- Controller can be attached left or right side
- Variable speed setting
- · Rear wheel drive
- Front wheel suspension



- Electromagnetic brake
- · Aeroplane safe battery
- Rechargeable Lithium-ion battery
- Anti-tip wheels
- Calf strap
- Lap belt
- Horn

SPECIFICATIONS

- Overall (L to unfolded footrest): 104.5 cm (L) x 55.5 cm (W) x 88.5 cm (H)
- Folded (L to folded footrest): 82 cm (L) x 27.5 cm (W) x 80 cm (H)
- Seat (H from floor): 42 cm (L) x 46 cm (W) x 49 cm (H)
- Backrest (H from seat): 44 cm (W) x 42 cm (H)
- Push Handle Height: 88 cm
- Turning Circle: 108 cm
- Speed: 6 km/hr
- Max Range: Up to 18 km
- Safe Climb Angle: 6°
- Max Climb Angle: 6°
- Motor: 2 x 250W high torque brushless motors
- Battery: 24V / 10.4Ah Lithium-Ion
- Charger: 2 Amp Lithium-Ion
- Tyres (flat free): Front: 16.5 cm (6.5") / Rear: 28 cm (11")
- Product Weight: 13.3 kg with battery, 12.3 kg without
- Max User Weight: 140 kg (on flat ground)

OPTIONAL ACCESSORIES

- 950007: IKon Carbon Fibre Controller Holder \$97.85
- 950023: Ikon Endura Attendant Joystick Mount for IK2PC or IK2PF \$35.00
- X-LUSB: Ikon Endura Multi-functional Light + Charging Accessory (LED Lamp, USB & Thunderport) - \$60.00
- 950016: Ikon Endura BLACK Overlay Kit (includes Back/Seat Cover & Leg Strap) \$95.00
- 950017: Ikon Endura BLUE Overlay Kit (includes Back/Seat Cover & Leg Strap) \$95.00
- 950018: Ikon Endura PURPLE Overlay Kit (includes Back/Seat Cover & Leg Strap) \$95.00
- 950019: Ikon Endura RED Overlay Kit (includes Back/Seat Cover & Leg Strap) \$95.00