



IKON ENDURO POWER LITE CARBON FIBRE ELECTRIC WHEELCHAIR

\$4,429.00



SKU: IK2PCS

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The ENDURO PowerLite Electric Wheelchair combines durability and convenience with its carbon fiber frame and foldable design for easy transport and storage. It features a comfortable padded seat and backrest, with a fold-down backrest that includes a handy storage pocket. Lift-up armrests and an onboard controller with an LCD display and fingertip controls ensure ease of use. The wheelchair boasts a battery power indicator, variable speed settings, and rear-wheel drive for smooth operation.

Enjoy a comfortable ride with front-wheel suspension, an electromagnetic brake, and anti-tip wheels. The airplane-safe, rechargeable Lithium-ion battery provides reliable power. Additional features include a calf strap, lap belt, and horn for added safety and convenience.

FEATURES

- Carbon fibre frame
- Folds up for transport and storage
- Comfortable padded seat and backrest
- Fold down backrest with storage pocket
- Lift up armrests
- Onboard controller with LCD display and fingertip controls- Battery power indicator
- Controller can be attached left or right side
- Variable speed setting
- Rear wheel drive
- Front wheel suspension
- Electromagnetic brake
- Aeroplane safe battery

- Rechargeable Lithium-ion battery
- Anti-tip wheels
- Calf strap
- Lap belt
- Horn

SPECIFICATIONS

- Overall (L to unfolded footrest): 104.5cm (L) x 55.5cm (W) x 88.5cm (H)
- Folded (L to folded footrest): 820cm (L) x 27.5cm (w) x 80cm (H)
- Seat (H from floor): 42cm (L) x 46cm (W) x 49cm (H)
- Backrest (H from seat): 44cm (W) x 42cm (H)
- Push Handle Height: 88cm
- Turning Circle: 108cm
- Speed: 6km/hr
- Max Range: Up to 18km²
- Safe Climb Angle: 6°
- Max Climb Angle: 6°
- Motor: 2 x 250W high torque brushless motors
- Battery: 24V / 10.4Ah Lithium-Ion
- Charger: 2 Amp Lithium-Ion
- Tyres (flat free): Front: 16.5cm (6.5") / Rear: 28cm (11")
- Product Weight: 13.3kg with battery, 12.3kg without
- Max User Weight: 140kg (on flat ground)