



IKON ENDURO POWER LITE FRONT FOLDING CARBON FIBRE POWERCHAIR - BLACK

\$4,429.00 Original price was:

\$4,429.00.~~**\$3,580.00**~~Current price is: **\$3,580.00**.Sale

Available Online Only!



SKU: IK2PCF

PRODUCT DESCRIPTION

Purpose & Benefits

The **Ikon Enduro Powerlite - Front Folding Carbon Fibre Electric Wheelchair** delivers the same trusted performance as Australia's best-selling carbon fibre wheelchair, now in a front-folding design. Built for travel, convenience, and independence, this model combines ultra-lightweight construction with powerful dual motors and exceptional comfort for everyday mobility.

Features

- Front-folding design for quick and compact transport and storage
- Lightweight carbon fibre frame – just 14.5 kg (without battery and cushions)
- Dual 250W high-torque brushless motors for smooth and reliable performance
- Airline-safe lithium-ion battery for easy travel
- Rear-wheel drive and front-wheel suspension for enhanced stability and comfort
- Onboard joystick controller with LCD display – attachable to either armrest
- Remote control included for flexible operation

- Padded seat and backrest for all-day comfort
- Folding footplate, anti-tip wheels, and rear storage pocket for added safety and convenience
- Supports users up to 150 kg

Specifications

Specification	Details
Overall Dimensions	900 mm (L) x 630 mm (W) x 970 mm (H)
Folded Dimensions	630 mm (L) x 340 mm (W) x 720 mm (H)
Seat Size	400 mm (W) x 400 mm (D) x 470 mm (H from floor)
Backrest Height	500 mm (from seat)
Push Handle Height	970 mm
Turning Radius	750 mm
Speed	Up to 6 km/h
Maximum Range	Up to 18 km*
Maximum Climb Angle	13°
Motor	2 x 250W High Torque Brushless Motors
Battery	24V / 10.4Ah Rechargeable Lithium-ion (Airline-safe)
Battery Weight	1.5 kg
Product Weight	14.5 kg (without battery and cushions)
Maximum User Weight	150 kg
Tyres	Flat-free – Front: 165 mm (6.5”) / Rear: 200 mm (7.9”)
Charger	2.0 Amp Lithium-Ion

Warranty

- 5-year manufacturer's warranty*

Included Items

- 1 x Ikon Enduro Powerlite Front Folding Electric Wheelchair
- 1 x 24V Lithium-ion Battery
- 1 x Battery Charger

- 1 x Remote Control
- 1 x User Manual
- 1 x Tool Kit

Note: Battery range varies based on user weight, terrain, and speed. The product is supplied fully assembled and ready to use.