

## IMMEDIA 4WAYGLIDE ADJUSTMENT MATTRESS



**SKU:** IM140/150



### PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

4WayGlide System consists of two parts that are fixed to the mattress and a Glide Mattress that lies on top. 4WayGlide in-bed system is suitable for heavy users with partial mobility in the legs.

4WayGlide Midi has a locking system between the nylon sheet and glide mattress and non-slip material along the sides of the nylon sheet.

### FEATURES

- Complete patient transfer bed system
- Ideal for those with reduced mobility
- Facilitates independent and assisted patient movement in bed
- Ideal for hospitals and care facilities
- Non-slip edges along sheet and mattress prevent patient slippage
- Mattress: Nylon back is low-static and low-friction for minimum pressure
- Maximum temperature for cleaning: 80°C
- Maximum temperature for tumble dry: 60°C
- Maximum user weight: 200 kg

### SPECIFICATIONS

- Top: Polyester/cotton
- Padding: Polyester fibre
- Underside: Nylon
- Handles: Polyester

