

IMMEDIA 4WAYGLIDE ADJUSTMENT MATTRESS



SKU: IM140/150



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

4WayGlide System consists of two parts that are fixed to the mattress and a Glide Mattress that lies on top. 4WayGlide in-bed system is suitable for heavy users with partial mobility in the legs.

4WayGlide Midi has a locking system between the nylon sheet and glide mattress and non-slip material along the sides of the nylon sheet.

FEATURES

- Complete patient transfer bed system
- Ideal for those with reduced mobility
- Facilitates independent and assisted patient movement in bed
- Ideal for hospitals and care facilities
- Non-slip edges along sheet and mattress prevent patient slippage
- Mattress: Nylon back is low-static and low-friction for minimum pressure
- Maximum temperature for cleaning: 80°C
- Maximum temperature for tumble dry: 60°C
- Maximum user weight: 200 kg

SPECIFICATIONS

- Top: Polyester/cotton
- Padding: Polyester fibre
- Underside: Nylon
- Handles: Polyester

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au