


## IMMEDIA EASYTURN SWIVEL CUSHION



**SKU:** IM990454

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Immedia EasyTurn Swivel Cushion is ideal for individuals with limited mobility. This clever, compact swivel cushion is designed to aid transfers to and from the bed, the car or a chair. The soft, flat cushion sits under the buttocks and helps the user to turn without strain to the back, hips or knees. The underside of the cushion is covered in a non-slip material to ensure the cushion doesn't slip while in use.

### FEATURES

- **Ease of Movement:** This cushion reduces strain on the knees, hips, and back by allowing a smooth, swivel motion that supports the user in pivoting to the desired position.
- **Non-Slip Base:** Equipped with a grip material on the underside, it firmly adheres to the seat, ensuring safety and stability during use.
- **Low-Friction Top Surface:** The top of the cushion is designed to be slick to decrease resistance, aiding in easier turning with minimal effort.
- **Portable and Lightweight:** Its compact design and light construction make it easy to transport and use in various settings – at home, in the car, or while travelling.
- **Supportive:** Provides a sturdy and comfortable platform for sitting, capable of accommodating a range of body types and weights.
- **Durable Material:** Made with quality materials to withstand regular use and maintain its effectiveness over time.

The Immedia EasyTurn Swivel Cushion is a practical solution for anyone who needs help turning when seated. Its simple design belies the profound impact it can have on the independence and comfort of individuals with mobility challenges, making daily tasks easier and safer.

## SPECIFICATIONS

Item	Product Code	Diameter (mm)	Material	Maximum user weight (kg)
EasyTurn Swivel	IM990454	450	Velour, polyester, nylon, polyurethane coating, polyethylene	150

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)