


## IMMEDIA MULTIGLIDE NYLON SLIDE SHEET WITH HANDLES



**SKU:** IM130DR-H

### PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Immedia MultiGlide Slide Sheet with handles is a tubular slide sheet with a high friction surface on the outside and an ultra-low friction material on the inside. The MultiGlide slide sheet is used during transfers to reduce friction at pressure points, minimising the risk of injury and shear. The waterproof slide sheet is ideal to facilitate transfers from a bed to an x-ray table or MT scanner and can also be used as a means for positioning, turning and moving patients up in bed. It is easy to apply and should always be removed after use.

Available in two sizes: Medium and Large.

### FEATURES

- Tubular high-low friction slide sheet with integrated handles for repositioning in bed
- Suitable for horizontal transfers
- Two strategically placed handles for user grip and handling
- High friction material on the outside and ultra-low friction material inside
- Made from durable waterproof Nylon
- Loop strap for easy storage and colour codes to differentiate between sizes
- Machine washable
- Maximum temperature for cleaning: 80 °C
- Maximum temperature for tumble dry: 60 °C
- Iron at a low temperature
- Maximum user weight: 300 kg

## SPECIFICATIONS

Item	Product Code	Width (mm)	Length (mm)	Material	Maximum user weight (kg)
MultiGlide w/handles - Medium (pink strap)	IM130DR	1300	800	Nylon with silicone and polyurethane coating, polyester	300
MultiGlide w/handles - Large (light blue strap)	IM160DR	1600	800	Nylon with silicone and polyurethane coating, polyester	300

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)