

IMMEDIA MULTIGLIDE NYLON SLIDE SHEET WITH HANDLES

\$269.00 - \$279.00 Price range: \$269.00 through \$279.00



SKU: IM130DR-H



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Immedia MultiGlide Slide Sheet with handles is a tubular slide sheet with a high friction surface on the outside and an ultra-low friction material on the inside. The MultiGlide slide sheet is used during transfers to reduce friction at pressure points, minimising the risk of injury and shear. The waterproof slide sheet is ideal to facilitate transfers from a bed to an x-ray table or MT scanner and can also be used as a means for positioning, turning and moving patients up in bed. It is easy to apply and should always be removed after use.

Available in two sizes: Medium and Large.

FEATURES

- Tubular high-low friction slide sheet with integrated handles for repositioning in bed
- Suitable for horizontal transfers
- Two strategically placed handles for user grip and handling
- High friction material on the outside and ultra-low friction material inside
- Made from durable waterproof Nylon
- Loop strap for easy storage and colour codes to differentiate between sizes
- Machine washable
- Maximum temperature for cleaning: 80 °C
- Maximum temperature for tumble dry: 60 °C
- Iron at a low temperature
- Maximum user weight: 300 kg

SPECIFICATIONS

Item	Product Code	Width (mm)	Length (mm)	Material	Maximum user weight (kg)
MultiGlide w/handles - Medium (pink strap)	IM130DR	1300	800	Nylon with silicone and polyurethane coating, polyester	300
MultiGlide w/handles - Large (light blue strap)	IM160DR	1600	800	Nylon with silicone and polyurethane coating, polyester	300