




IMMEDIA TRANSFER BOARD

\$135.00 - \$1,000.00 Price range: \$135.00 through \$1,000.00



SKU: IM40-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Immedia Butterfly Transfer Board is a long, wide and stable transfer board with a smooth top surface and special cutouts for wheelchair positioning. The fibreglass Butterfly board is ideal for bridging small gaps during seated transfers. Equipped with non-slip pads on the underside that prevent the board from sliding out of place, the board allows for good positioning during transfers. With a larger sitting and sliding area, the Butterfly board is helpful for transfers made in short advances.

FEATURES

- A long, wide and stable transfer board
- Bridges short gaps during seated transfers
- Special cutouts for wheelchair positioning
- Non-slip pads on the underside for stability
- GlideWing facilitates smooth seated transfers (additional slide function)
- Available with large or extra large GlideWing sold separately
- Wipe down with an alcohol wipe
- Maximum user weight: 250 kg

SPECIFICATIONS

| Item | Product Code | Length (mm) | Width (mm) | Material | Product weight (kg) | Maximum user weight (kg) |
|------|--------------|-------------|------------|----------|---------------------|--------------------------|
|------|--------------|-------------|------------|----------|---------------------|--------------------------|

| | | | | | | |
|----------------------------------|--------|-----|-----|---------------------|-----|-----|
| Butterfly Board | IM401 | 660 | 320 | Fibreglass | 2,2 | 250 |
| Butterfly GlideWing L | IM4001 | 400 | 350 | Nylon, polyester | 2,2 | 250 |