




IMMEDIA TRANSFER BOARD

\$135.00 - \$1,000.00 Price range: \$135.00 through \$1,000.00



SKU: IM40-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Immedia Butterfly Transfer Board is a long, wide and stable transfer board with a smooth top surface and special cutouts for wheelchair positioning. The fibreglass Butterfly board is ideal for bridging small gaps during seated transfers. Equipped with non-slip pads on the underside that prevent the board from sliding out of place, the board allows for good positioning during transfers. With a larger sitting and sliding area, the Butterfly board is helpful for transfers made in short advances.

FEATURES

- A long, wide and stable transfer board
- Bridges short gaps during seated transfers
- Special cutouts for wheelchair positioning
- Non-slip pads on the underside for stability
- GlideWing facilitates smooth seated transfers (additional slide function)
- Available with large or extra large GlideWing sold separately
- Wipe down with an alcohol wipe
- Maximum user weight: 250 kg

SPECIFICATIONS

Item	Product Code	Length (mm)	Width (mm)	Material	Product weight (kg)	Maximum user weight (kg)
------	--------------	-------------	------------	----------	---------------------	--------------------------

Butterfly Board	IM401	660	320	Fibreglass	2,2	250
Butterfly GlideWing L	IM4001	400	350	Nylon, polyester	2,2	250