


IMMEDIA TRANSFER BOARD



SKU: IM40-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Immedia Butterfly Transfer Board is a long, wide and stable transfer board with a smooth top surface and special cutouts for wheelchair positioning. The fibreglass Butterfly board is ideal for bridging small gaps during seated transfers. Equipped with non-slip pads on the underside that prevent the board from sliding out of place, the board allows for good positioning during transfers. With a larger sitting and sliding area, the Butterfly board is helpful for transfers made in short advances.

FEATURES

- A long, wide and stable transfer board
- Bridges short gaps during seated transfers
- Special cutouts for wheelchair positioning
- Non-slip pads on the underside for stability
- GlideWing facilitates smooth seated transfers (additional slide function)
- Available with large or extra large GlideWing sold separately
- Wipe down with an alcohol wipe
- Maximum user weight: 250 kg

SPECIFICATIONS

Item	Product Code	Length (mm)	Width (mm)	Material	Product weight (kg)	Maximum user weight (kg)
------	--------------	-------------	------------	----------	---------------------	--------------------------

Butterfly Board	IM401	660	320	Fibreglass	2,2	250
Butterfly GlideWing L	IM4001	400	350	Nylon, polyester	2,2	250

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au