



INDEELIFT PPU-S - PEOPLE PICKER UPPER STANDING LIFT



SKU: HFL-PPU-S

PRODUCT DESCRIPTION

PURPOSE & BENEFITS OF INDEELIFT PEOPLE PICKER UPPER STANDING LIFT

The Indeelift People Picker Upper Standing Lift is a groundbreaking mobility aid designed to safely and efficiently lift individuals from the floor to a standing or seated position. This innovative device addresses a critical need in both home and healthcare settings, offering a dignified solution for those who have fallen or need assistance getting up. It's particularly beneficial for the elderly, people with disabilities, or those recovering from surgery, reducing the risk of injury for both the user and caregivers.

FEATURES

- **Safe and Secure Lifts:** The Indeelift is engineered to provide a stable and secure lift, minimizing the risk of falls or further injury when attempting to stand or transfer to a wheelchair or other seating.
- **User-Operated:** Designed to allow users to operate the lift independently if able, promoting self-reliance and confidence. This feature is particularly empowering for individuals who live alone or wish to reduce dependence on others for assistance with getting up.
- **Adjustable Height:** The device can adjust to various height settings, accommodating different user heights and preferences, and facilitating transfers to a wide range of seating surfaces or standing positions.
- **Compact and Portable Design:** While robust, the Indeelift has a compact footprint and is designed for easy maneuverability within homes or healthcare facilities. Some models may also be portable, offering the flexibility to use the device in different settings.
- **Heavy-Duty Construction:** Built with durable materials, the Indeelift is capable of lifting individuals up to a certain weight limit, making it a suitable option for a diverse range of users,

including bariatric individuals.

- **Simple to Use:** The lift's controls are intuitive, ensuring ease of use for individuals or caregivers, without the need for extensive training.

The Indeelift People Picker Upper Standing Lift represents a significant advancement in mobility and safety for individuals at risk of falls or those needing assistance to stand. By offering an innovative and humane solution to a common challenge, the Indeelift improves the quality of life, independence, and safety for users and caregivers alike.

SPECIFICATIONS

- Safe Working Load: 180kg
- Overall Depth: 26.5 " (67 cm)
- Maximum Overall Height (raised): 70.75" (174 cm)
- Minimum Overall Height (lowered): 41.75" (100 cm)
- External Width incl Rise Handles: 23" (58 cm)
- Maximum Seat Height: 30" (76 cm)
- Seat Width: 14" (35 cm)
- Wheel Diameter: 5" (13 cm)
- Unit Weight: 78 lbs. (35 kg)

BROCHURE

[PPU-S - Data Sheet](#)

VIDEO

[embed]<https://www.youtube.com/watch?v=58YwEKMhS5k>[/embed]

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au